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Bruker D8 Advance Technical Specification

This instrument is based around a two circle goniometer, enclosed in a radiation safety enclosure.

X-Ray Source:

The X-ray Source is a 2.2kW Cu anode long fine focus ceramic X-ray tube. (Mo, Co and Cr also available). The running conditions for the X-Ray tube are 40 kV and 40 mA, the power supply controls these with a stability of better than 0.01% for the high voltage and 10% of the variation of the supply for the current.

Goniometer

This is a high precision microprocessor controlled, two circle goniometer with independent stepper motors and optical encoders for the Theta and 2Theta circles with the following specifications:

Maximum measurement circle diameter
250 mm

Minimum measurement circle diameter
100 mm

Smallest angular step size 0.0001°

Reproducibility +/- 0.0001°

Maximum rotational speed 1500°/min

Angular range (Theta) -5° to 40°

Angular range (2Theta) -10° to 60°

Detector

The detector is a NaI dynamic scintillation detector with a maximum count rate of $2 \times 10^6 \text{ s}^{-1}$ (although it should not be exposed to in excess of $5 \times 10^5 \text{ s}^{-1}$ for periods longer than about 1 second). A computer controlled absorber is mounted directly in front of the detector and is used to attenuate the beam by about 2 orders of magnitude. To change the angular resolution several slits (0.05 to 2 mm) are available to be placed in front of the detector, as are a set of Soller slits to allow for small samples. Attenuators made of Cu and Ni are also available and can be manually placed in the beam, by inserting them into the slit holders in front of the detector. The detector electronics are capable of a count rate up to $3 \times 10^7 \text{ s}^{-1}$.

Software

The software suit is all windows based, and enables the user to both acquire and analyse data. Data acquisition is possible both directly by the immediate measurement program "Adjust" and the automated "Job" program via a DQL file. In the later case the data can be displayed in real time using the display program.

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